



<b>Period of Placement in Situ</b>	Short and medium duration (1-3 days)	Short and medium duration (1-3 days)	Long duration (up to 7 days)																																																																
<b>Typical Clinical Application</b>	Flushing obstructed uroliths retrograde into the bladder	Securing urine passage. Typical used in FLUDD/FIC cats	Hospitalized patients and patients undergoing surgery																																																																
<b>Features and Benefits</b>	<p>The cat catheter with open end allows the vet or technician to flush uroliths retrograde into the bladder</p> <p>The catheter is more rigid, than the catheters with side holes</p> <p>Suturing holes in the baseplate allows suturing to the prepuce, but will easier get obstructed than the catheters with wide holes, as this catheter only have one hole</p> <p>Luer lock fitting: for easy connection to syringe or urine collection bag</p>	<p>The catheter has two lateral side holes and rounded, atraumatically closed end</p> <p>Comes in different sizes and with and without stylet</p> <p>Suturing holes in the baseplate allows suturing to the prepuce</p> <p>Luer lock fitting: for easy connection to syringe or urine collection bag</p>	<p>BUSTER Barium catheter is made of polyurethane and contains Barium</p> <p>The Barium makes it visible in x-ray</p> <p>Polyurethane softens at body temperature, which increases the biocompatibility, hence suitable for longer term placement and hospitalized patients</p> <p>The catheter has two lateral side holes and rounded, atraumatically closed end</p> <p>Suturing holes in the baseplate allows suturing to the prepuce</p> <p>Luer lock fitting: for easy connection to syringe or urine collection bag</p>																																																																
<b>Catheter Types</b>	<b>BUSTER Cat Catheter, open end</b>	<b>BUSTER Cat Catheter, side holes</b>	<b>BUSTER Barium Cat Catheter</b>																																																																
<b>Sizes and Packaging</b>	<table border="1"> <thead> <tr> <th>Cat. No</th> <th>Description</th> <th>Size</th> <th>Package</th> </tr> </thead> <tbody> <tr> <td>273338</td> <td>Without stylet</td> <td>1.0 x 130 mm / 3 Fr x 5½ in</td> <td>12/pk</td> </tr> <tr> <td>273339</td> <td>Without stylet</td> <td>1.3 x 130 mm / 4 Fr x 5½ in</td> <td>12/pk</td> </tr> </tbody> </table> <p><b>BUSTER Easy Slide Cat Catheter, open end</b></p> <table border="1"> <tbody> <tr> <td>273436</td> <td>Without stylet</td> <td>1.2 x 110 mm / 3.5 Fr x 4½ in</td> <td>5/pk</td> </tr> <tr> <td>273437</td> <td>Without stylet</td> <td>1.2 x 140 mm / 3.5 Fr x 5½ in</td> <td>5/pk</td> </tr> <tr> <td>273438</td> <td>Without stylet</td> <td>1.2 x 180 mm / 3.5 Fr x 7 in</td> <td>5/pk</td> </tr> </tbody> </table>	Cat. No	Description	Size	Package	273338	Without stylet	1.0 x 130 mm / 3 Fr x 5½ in	12/pk	273339	Without stylet	1.3 x 130 mm / 4 Fr x 5½ in	12/pk	273436	Without stylet	1.2 x 110 mm / 3.5 Fr x 4½ in	5/pk	273437	Without stylet	1.2 x 140 mm / 3.5 Fr x 5½ in	5/pk	273438	Without stylet	1.2 x 180 mm / 3.5 Fr x 7 in	5/pk	<table border="1"> <thead> <tr> <th>Cat. No</th> <th>Description</th> <th>Size</th> <th>Package</th> </tr> </thead> <tbody> <tr> <td>273461</td> <td>Without stylet</td> <td>1.3 x 130 mm / 4 Fr x 5½ in</td> <td>12/pk</td> </tr> <tr> <td>273453</td> <td>With stylet</td> <td>1.0 x 130 mm / 3 Fr x 5½ in</td> <td>12/pk</td> </tr> <tr> <td>273463</td> <td>With stylet</td> <td>1.3 x 130 mm / 4 Fr x 5½ in</td> <td>12/pk</td> </tr> <tr> <td>273456</td> <td>With stylet</td> <td>1.3 x 160 mm / 4 Fr x 6½ in</td> <td>12/pk</td> </tr> </tbody> </table> <p><b>BUSTER Easy Slide Cat Catheter, side holes</b></p> <table border="1"> <tbody> <tr> <td>273400</td> <td>Without stylet</td> <td>1.2 x 110 mm / 3.5 Fr x 4½ in</td> <td>5/pk</td> </tr> <tr> <td>273401</td> <td>Without stylet</td> <td>1.2 x 140 mm / 3.5 Fr x 5½ in</td> <td>5/pk</td> </tr> </tbody> </table>	Cat. No	Description	Size	Package	273461	Without stylet	1.3 x 130 mm / 4 Fr x 5½ in	12/pk	273453	With stylet	1.0 x 130 mm / 3 Fr x 5½ in	12/pk	273463	With stylet	1.3 x 130 mm / 4 Fr x 5½ in	12/pk	273456	With stylet	1.3 x 160 mm / 4 Fr x 6½ in	12/pk	273400	Without stylet	1.2 x 110 mm / 3.5 Fr x 4½ in	5/pk	273401	Without stylet	1.2 x 140 mm / 3.5 Fr x 5½ in	5/pk	<table border="1"> <thead> <tr> <th>Cat. No</th> <th>Description</th> <th>Size</th> <th>Package</th> </tr> </thead> <tbody> <tr> <td>273467</td> <td>With stylet</td> <td>1.0 x 130 mm / 3 Fr x 5½ in</td> <td>12/pk</td> </tr> <tr> <td>273468</td> <td>With stylet</td> <td>1.3 x 130 mm / 4 Fr x 5½ in</td> <td>12/pk</td> </tr> </tbody> </table>	Cat. No	Description	Size	Package	273467	With stylet	1.0 x 130 mm / 3 Fr x 5½ in	12/pk	273468	With stylet	1.3 x 130 mm / 4 Fr x 5½ in	12/pk
Cat. No	Description	Size	Package																																																																
273338	Without stylet	1.0 x 130 mm / 3 Fr x 5½ in	12/pk																																																																
273339	Without stylet	1.3 x 130 mm / 4 Fr x 5½ in	12/pk																																																																
273436	Without stylet	1.2 x 110 mm / 3.5 Fr x 4½ in	5/pk																																																																
273437	Without stylet	1.2 x 140 mm / 3.5 Fr x 5½ in	5/pk																																																																
273438	Without stylet	1.2 x 180 mm / 3.5 Fr x 7 in	5/pk																																																																
Cat. No	Description	Size	Package																																																																
273461	Without stylet	1.3 x 130 mm / 4 Fr x 5½ in	12/pk																																																																
273453	With stylet	1.0 x 130 mm / 3 Fr x 5½ in	12/pk																																																																
273463	With stylet	1.3 x 130 mm / 4 Fr x 5½ in	12/pk																																																																
273456	With stylet	1.3 x 160 mm / 4 Fr x 6½ in	12/pk																																																																
273400	Without stylet	1.2 x 110 mm / 3.5 Fr x 4½ in	5/pk																																																																
273401	Without stylet	1.2 x 140 mm / 3.5 Fr x 5½ in	5/pk																																																																
Cat. No	Description	Size	Package																																																																
273467	With stylet	1.0 x 130 mm / 3 Fr x 5½ in	12/pk																																																																
273468	With stylet	1.3 x 130 mm / 4 Fr x 5½ in	12/pk																																																																
<b>Related Products</b>	KRUUSE Urine strips, Cat. No 290798, 100/pk	KRUUSE Urine strips, Cat. No 290798, 100/pk BUSTER Urine bags, Luer Lock, Cat. No 273855, 10/pk BUSTER Hanger for Urine Bag, Cat. No 273857, 10/pk	KRUUSE Urine strips, Cat. No 290798, 100/pk BUSTER Urine bags, Luer Lock, Cat. No 273855, 10/pk BUSTER Hanger for Urine Bag, Cat. No 273857, 10/pk																																																																