

Comfortable, color coded muzzles

BUSTER Easy-ID Dog Nylon Muzzle

- Extra padding around the nose
- Color coded for easy size identification
- Reinforced stitching to ensure long-lasting use
- Available in 10 sizes from 2XS to 5XL



Selecting the right muzzle for the dog

- Use a string to measure around the dog's nose at the position shown.
- Place your index finger under the string to ensure enough space for the dog to be comfortable. Do not tighten the string; it should be free to move around the dog's nose.
- Lay the string along a ruler to get appropriate length, then select the proper muzzle. The muzzle must not be tight. It is better a little loose than too tight.



Weight/Breed Guideline

Covetrus #	KRUUSE #	Size	Patient Weight Range (lbs)	Snout Circumference (")	Colour	Breeds
068205	279465	2XS	< 5 lbs	9 cm / 3.5"	●	Chihuahua
068206	279466	XS	5-10 lbs	11.5 cm / 4.5"	●	Toy Poodle
068194	279467	S	10-25 lbs	13 cm / 5"	●	Dachshound
068195	279468	M	25-40 lbs	15 cm / 6"	●	Cocker Spaniel
068196	279469	L	40-70 lbs	16.5 cm / 6.5"	●	Bull Terrier
068197	279470	XL	70-100 lbs	20 cm / 8"	●	Doberman
068198	279471	2XL	100-150 lbs	25 cm / 10"	●	German Shepherd
068199	279472	3XL	100-120 lbs	28 cm / 11"	●	Labrador Retriever
068200	279473	4XL	120-150 lbs	30 cm / 12"	●	Great Dane
068201	279475	5XL	> 150 lbs	35 cm / 14"	●	Giant Breeds
068202	279477	XXS-5XL	Set of 10 BUSTER Easy-ID Dog Nylon Muzzles		Ass.	All
068204	279487	S-2XL	Set of 5 BUSTER Easy-ID Dog Nylon Muzzles		Ass.	All
068210	279489	XS - 2XL	Set of 5 BUSTER Easy-ID Dog Nylon Muzzles + 1 BUSTER muzzle for cat (size M)		Ass.	All